

ROHO®

# Hybrid Select



permobil

# Exceptional seating with hybrid positioning

The Hybrid Select pelvic loading area adjusts to your unique shape, weight, and positioning changes. The ROHO® overlay and IT insert utilizes the ISOFLO® Memory Control which allows to find and set the right position each and every time.

## Unique lateral stability

The innovative design of the Hybrid Select uses the positioning capability of a contoured foam base with high resiliency (HR) foam and leg troughs to offer you support, lateral stability, and provide relief to the coccyx/tailbone.

## Proven skin protection

With Hybrid Select, experience truly functional stability combined with great skin and soft tissue protection. The contoured design of the cushion allows the pelvis to be immersed in air which shapes and moves to the body, and protects the skin and soft tissue.

## Optional offloading comes standard

The Hybrid Select provides flexibility and forgiveness that other offloading cushions do not, making it easier to be in the right position. It offers an exclusive optional offloading technique by removing the ischial tuberosities (IT) air insert.

## Attributes of great seating



### Skin Protection

Immersion test (ISO 16840-2) reveals how deeply the body sinks into the cushion.

Envelopment test (ISO 16840-12) reveals how the body weight is cradled by the cushion, using envelopment or offloading.

Shear test (Shear force sensor) reveals the potential of the cushion interaction with the body to distort/deform the tissue.

Contact Area (Pressure mapping) reveals how much contact is being made between the seat support and the individual.



### Positioning

Slide resistant test (ISO 16840-2) reveals how much force is required to slide forward on the cushion when seated.




### Stability

Lateral tilt test (ISO 16840-13) reveals the stability the cushion may provide during side-to-side leaning.



Good posture with healthy weight distribution for comfortable seating



Exceptional relief and comfort


# Hybrid Select Specifications

Cover type	Stretch-Air®
Product weight	1,45 kg*
Max. user weight	227 kg
Cushion warranty	36 month limited
Cover warranty	6 month limited



## Size options

35.5 cm x 37.5 cm	40.5 cm x 42.5 cm	43.0 cm x 47.5 cm	51.0 cm x 42.5 cm
35.5 cm x 42.5 cm	40.5 cm x 45.0 cm	43.0 cm x 52.5 cm	51.0 cm x 45.0 cm
38.0 cm x 40.0 cm	40.5 cm x 47.5 cm	45.5 cm x 42.5 cm	51.0 cm x 47.5 cm
38.0 cm x 42.5 cm	40.5 cm x 52.5 cm	45.5 cm x 45.0 cm	51.0 cm x 52.5 cm
38.0 cm x 45.0 cm	43.0 cm x 42.5 cm	45.5 cm x 47.5 cm	
40.5 cm x 40.0 cm	43.0 cm x 45.0 cm	45.5 cm x 52.5 cm	




Hybrid Select  
Front view



Hybrid Select  
left side view



Hybrid Select  
back view



\*Product weight is based on a size 40,5 x 40,5 cushion. Product weight varies depending on cushion size. Measurements are approximates and subject to change without notice.

The Hybrid Select gives me the positioning I need and makes it easier for me to transfer.

-Eric, Permobil user

ROHO®

## Hybrid Select

The Hybrid Select is an advanced cushion providing exceptional positioning and stability, in addition to the ability to offload.

It uses the positioning capabilities of a contoured, high resiliency foam base with deep leg troughs to offer support, lateral stability, and provide relief to the coccyx.

**permobil**

permobil.com | info@permobil.co.uk | T. +44 (0)1484 722 888



HYBRIDSELECTDATASHEETUK211206